

BEAUTIFIL II LS

CLICK / SCAN



FOR MORE INFO

UNIVERSAL BIOACTIVE NANO-HYBRID LOW SHRINKAGE COMPOSITE RESIN

Designed for clinicians who prefer a stiffer, non-sticky, sculptable composite that offers...

- **Predictable Aesthetics and Function** to meet patient expectations
- **Easy to Sculpt, Carve & Polish**
- **Eliminates Post-op Sensitivity &** protect against secondary caries

Formulated with ground breaking SRS Monomer Technology & patented S-PRG filler, **BEAUTIFIL II LS** exhibits extremely low polymerization shrinkage with bioactive properties to ensure long lasting aesthetic restorations.

Simplified shade line-up helps minimize inventory and achieve predictable aesthetics with simple or complex composite layering techniques.



Available in 14 Shades – 4gm syringe

- Universal Dentin – A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3
- Opaque Dentin – A20, A30
- Enamel – INC, Bleach White - BW

CLINICAL APPLICATIONS

- Direct veneers
- Cosmetic tooth modifications & diastema closure
- Class I – VI cavities and root caries
- Tooth wear – abrasion, erosion, attrition, abfraction
- Repair of direct and indirect restorations

FEATURES

- ✓ **Aesthetics that mimic natural teeth**
 - Easy shade match with 3 different opacities & opalescence
 - Shade stability before and after cure
 - Effortless polish to natural lustre
- ✓ **Bioactive to protect against secondary caries**
 - Minimize plaque formation
 - Reduce cariogenic bacteria
- ✓ **Superior handling**
 - Non-sticky, non-tacky & sculptable
 - Easy to adapt and place
- ✓ **Long lasting restorations**
 - Low polymerization shrinkage
 - Strong and durable
 - Remarkable stain & wear resistance

PROPERTIES

- Compressive Strength: 364 Mpa
- Polymerization Shrinkage: 0.85% (vol)
- Filler Load: 83 wt%

RELATED PRODUCTS

- Beautiful II Enamel & Gingiva
- BeautiBond / FL Bond II
- Super-Snap Xtreme Disks / OneGloss 2-in-1 Smart Polisher
- Diamond Polishing Paste with Super-Snap Buff



Official Partner

