

# BEAUTIFIL Injectable X

CLICK / SCAN



FOR MORE INFO

## ONE TRUE BIOACTIVE UNIVERSAL INJECTABLE RESTORATIVE

Designed to help clinicians simplify work flow and save time with predictability...

- **Inject as you Shape** to boost efficiency
- **Remarkable Strength** for anterior and posterior restorations
- **Life-like Aesthetics** with self-polishing gloss

Formulated with the Nano-optimised Silane technology to provide strength and aesthetics of a packable composite with the unique “inject as you shape” convenience.

Custom designed, disposable Barrier Sleeves for Beautiful Injectable X provides additional protection against cross-contamination.



### Available in 15 Shades – 2.2gm syringe

- Universal Dentin – SA0.5, A1, A2, A3, A3.5, A4, B1, B2, C2, D2
- Opaque Dentin – SA0.5O, A1O, A2O
- Enamel – INC, Bleach White - BW

## CLINICAL APPLICATIONS

- Restoration of Anterior / Posterior teeth...
  - Class I – VI cavities & Root caries
  - Pediatric and geriatric dentition
  - Tooth wear (Abrasion, Erosion, Abfraction)
- Minimally invasive cosmetic restorations
  - Diastema closure, Tooth surface defects
  - Tooth size and shape modification
  - Smile rejuvenation
- Cuspal Restoration
- Post-ortho corrections & Black triangle closure
- Composite Tooth Splinting
- Repair of Composite or Ceramic restorations

## FEATURES

- ✓ **Aesthetics that mimic natural teeth**
  - Easy shade match with 3 different opacities & opalescence
  - Seamless blend with natural tooth margins
  - Effortless gloss with a self-polish capability
- ✓ **Bioactive to protect against secondary caries**
  - Minimize plaque formation, reduce cariogenic bacteria & neutralize sugar acids
  - Fluoride release & recharge
- ✓ **Superior handling**
  - “Inject as you shape” convenience
  - Fast & efficient precise placement
- ✓ **Long lasting restorations**
  - Strong & durable with optimal filler-resin ratio
  - Remarkable stain & wear resistance

## PROPERTIES

- Compressive Strength: above 400 Mpa
- Flexural Strength: 127 Mpa
- Radiopaque

## RELATED PRODUCTS

- Beautiful Injectable X<sub>SL</sub>
- Barrier Sleeves
- BeautiBond / FL Bond II
- Super-Snap Xtreme Disks / OneGloss 2-in-1 Smart Polisher



Official Partner

