

BEAUTIFIL Injectable X^{SL}

CLICK / SCAN



FOR MORE INFO

ONE TRUE BIOACTIVE SELF-LEVELLING INJECTABLE RESTORATIVE

Designed to help clinicians simplify work flow and save time with predictability...

- **High strength self-levelling flowable**
- **Life-like Aesthetics** with self-polishing gloss
- **Durable** even on the surface of restorations

Formulated with the Nano-optimised Silane technology to provide strength and aesthetics of a packable composite with the unique flowable consistency.

Custom designed, disposable Barrier Sleeves for Beautiful Injectable X^{SL} provides additional protection against cross-contamination.



Flow after 60 sec

Available in 13 Shades – 2.2gm syringe

- Universal Dentin – A1, A2, A3, A3.5, A4, B1, B2
- Opaque Dentin – A2O, A3O
- Enamel – INC, Bleach White - BW
- Cervical – CV, Milky - MI

CLINICAL APPLICATIONS

- Restoration of Anterior / Posterior teeth...
 - Class III, IV, V cavities & Root caries
 - Gingival box of Class II
 - Tooth wear (Abrasion, Erosion, Abfraction)
 - Deep Margin Elevation (DME)
- Palatal shell technique for Class IV
- Base & Liner
- Orthodontic bracket placement
- Composite Tooth Splinting

FEATURES

- ✓ **Aesthetics that mimic natural teeth**
 - Easy shade match with 3 different opacities
 - Stable shades pre & post cure
 - Seamless blend with natural tooth margins
 - Effortless gloss with a self-polish capability
- ✓ **Bioactive to protect against secondary caries**
 - Minimize plaque formation, reduce cariogenic bacteria & neutralize sugar acids
 - Fluoride release & recharge
- ✓ **Superior handling**
 - Self-levelling flowable consistency
 - Fast & efficient precise placement
- ✓ **Long lasting restorations**
 - Strong and durable
 - Remarkable stain & wear resistance

PROPERTIES

- Compressive Strength: above 400 Mpa
- Flexural Strength: 127 Mpa
- Radiopaque

RELATED PRODUCTS

- Beautiful Injectable X, Beautiful II, Beautiful II LS
- Barrier Sleeves
- BeautiBond / FL Bond II
- Super-Snap Xtreme Disks / OneGloss 2-in-1 Smart Polisher



Official Partner

