

CLICK / SCAN



FOR MORE INFO

BEAUTIFIL OPAQUER

UNIVERSAL BIOACTIVE NANO-HYBRID FLOWABLE OPAQUER


Designed to effectively mask discoloration with a thin layer to help create invisible restorations...

- A **“Must Have”** for masking discolorations, stains, spots and even metal
- **Ideal Flowability** to achieve a uniform bubble-free thin layer

Formulated as a flowable opaque composite with the proven S-PRG filler technology to provide bioactive benefits while effectively masking deep discoloration with a thin layer. Reputed as the secret to success when creating minimally invasive life-like restorations.



Available in 2 Shades – 2gm syringe



CLINICAL APPLICATIONS

- Mask tooth discoloration due to
 - Internal staining
 - Amalgam tattoo
 - Trauma
 - Pulp capping
- Intraoral repair of PFM restorations
- Mask abutment tooth discoloration for indirect restorations

FEATURES

- ✓ **Optimal opacity to effectively mask discolorations**
 - 2 unique shades
 - Effective even with a very thin layer
- ✓ **Bioactive to protect against secondary caries**
 - Minimize plaque formation
 - Reduce cariogenic bacteria
- ✓ **Superior handling**
 - Ideal viscosity for uniform bubble-free, paint on application
- ✓ **Unique Syringe Design**
 - 360° multidirectional grip
 - Plunger designed to minimize wastage with pull-back mechanism

PROPERTIES

- Fluoride release and recharge
- Radiopaque

RELATED PRODUCTS

- Beautifil II / Beautifil II LS/ Beautifil Injectable X
- BeautiBond / FL Bond II
- Universal Primer



Official Partner

