

# BEAUTIFIL-Bulk

CLICK / SCAN



FOR MORE INFO

## BIOACTIVE BULK FILL RESTORATIVE


Designed to add convenience and predictability for time-efficient posterior restorations...

- **Faster Easier & Simpler Restorations** with a single increment of 4mm
- **Simplified Inventory** with just 2 shades that offer a remarkable shade match
- **Bioactive Benefits** to prevent secondary caries

Formulated with a unique balance of filler- resin chemistry with new monomers and a filler complex that includes the patented S-PRG filler technology to provide, a deep cure, minimal polymerization shrinkage stress and optimal physical properties to ensure lasting restorations that blend well with natural teeth.



Available in 2 Shades – 4.5gm syringe



## CLINICAL APPLICATIONS

- Restoration of Posterior teeth...
  - Class I and II cavities
  - Cusp build-up
- Core build-up
- Deciduous restorations

## FEATURES

- ✓ **Sculptable handling for bulk filling**
  - Deep cure in a single step upto 4mm
  - Fast & efficient enamel and dentin replacement
  - Non sticky and easy to adapt
- ✓ **Mimics the shade of the natural tooth**
  - 2 versatile shades that blends well
  - Stable shades pre & post cure
  - Ideal optical properties with remarkable shade match
- ✓ **Bioactive to protect against secondary caries**
  - Minimize plaque formation, reduce cariogenic bacteria & neutralize sugar acids
  - Fluoride release & recharge
- ✓ **Long lasting restorations**
  - Strong and durable
  - Minimal polymerization shrinkage
  - Increased wear resistance

## PROPERTIES

- Compressive Strength: above 340Mpa
- Flexural Strength: 113 Mpa
- Polymerisation shrinkage: 1.7%
- Filler load: 87 wt%
- Radiopaque

## RELATED PRODUCTS

- Beautiful Bulk Flow , Beautiful Injectable X, XSL
- BeautiBond / FL Bond II
- Dura White Stone, CompoMaster
- OneGloss 2-in-1 Smart Polisher



Official Partner

