



ALL PURPOSE TRUE-COLOUR STAIN SYSTEM

Designed for individual characterization of all types of ceramics to achieve life-like 2D or 3D stain effects...

- Low-fusing stains compatible with all conventional and CAD/CAM dental ceramics
- Ultra-fine glass technology assures easy, uniform application
- Intense colour pigment structure with natural fluorescence
- True-colour liquid to visualize results in real time

Formulated as an innovative, universal, vivid stain system that ensures predictable staining and characterization of all types of ceramics.



Available in sets and refills:

- Vintage Art Universal Basic Set
- Stain Refills – 2gm
- Glazing Powder GP – 15gm, 50gm
- Glazing Powder GP-F – 15gm, 50gm
- Vintage Art Universal Liquid, 50ml
- Yamamoto Liquid, 50ml

CLINICAL APPLICATIONS

- Staining of monolithic or layered ceramics
- Individualizing low & high fusing ceramics
- Ideal for both 2D & 3D staining
- Internal or external characterization

FEATURES

- ✓ **Convenience of Use**
 - Microfine particle structure that assures ease of mixing and application in thin, even and bubble-free layers to create different contrasts
- ✓ **Universal application – Streamlined Inventory**
 - Compatible with monolithic and layered/ High or Low-fusing ceramics
- ✓ **Flexibility**
 - Option of 2 glaze powders – with or without fluorescence
- ✓ **Time-efficient, Predictable & Error-free**
 - True-colour (Yamamoto) liquid with glass-like refractive index to visualize shade effect and translucency of the stain mixture before firing

FINAL FIRING TEMPERATURE GUIDELINE

- Veneered Lithium Disilicate
730°C – 740°C
- Monolithic Lithium Disilicate
740°C – 750°C
- Veneered Zirconia / PFM Restorations
750°C – 800°C
- Monolithic Zirconia / Porcelain Teeth
800°C – 850°C

RELATED PRODUCTS

- Vintage Pro
- Vintage ZR
- Vintage LD
- SHOFU ZR Lucent
- SHOFU Porcelain Brush Set